

## Accessories for shower/commode chair M2

### Cross bar and cross guard



#### Crossbar

↑ The crossbar can be placed on the two armrests to give extra support to the user. Is fastened with the two rubber bands. It gives the user feels like more safety.

**Cannot be used on M2, 200 kg's, 300 kg's and M2 mini.**



#### Cross guard

↑ The cross guard can be placed so it gives support in front of the user.

**Cannot be used on M2 200 kg and M2 300 kg's, M2 small and M2 mini.**



↑ Holder for side supports/Cross guard

Holders for cross guard are standard on M2-tip. For other chairs you order this mounting set, which is screwed on with the help of the screws for the back rest. Cross guard holders can be mounted together with holders for head support.



**HMN** since 1934

<b>Technical data</b>	
<b>Crossbar</b>	
Total length:	62,5 cm
Diameter of foam part:	App. 36 mm
Metal:	Stainless Steel
Rubber bands:	Rubber
Weight:	0,8 kg
Article no. (Cannot be used on M2, 200 kg's and M2 small)	310228
Article no. The Crossbar to M2 Small	310360
<b>Cross guard</b>	
Total width:	62,5 cm
Depth: (from back rest to centre of guard)	App. 37 cm
Diameter of foam part:	App.34 mm
Metal:	Stainless Steel
Weight :	1,6 kg
Article no. (Cannot be used on M2 200 kg's and M2 small)	310229
For ordinary M2's you need a set of holders article no. 310227 as well. On M2-Tip-Chairs holders on the chairs are standard. If you use a cross guard you cannot have side supports on the chairs as well.	
<b>Side support/Cross guard holders</b>	
Set with 2 pcs. which are mounted with the enclosed screws.	
Metal:	Stainless Steel
Weight (set with 2 pcs.):	0,5 kg
Holders are standard on M2-tip. For other chairs order as underneath	
Article no. (set with 2 pcs)	310227
Article no. One for right side	310230
Article no. One for left side	310231



**HMN** since 1934