

Accessories for shower/commode chair M2

Leg supports and heel stop



The leg support is made from net weave and can be adjusted by velcro around the foot rests. The leg support can be used on all M2 shower/commode chairs, and it exists in two sizes: standard and XL.

The leg support is recommended for all the tip-chairs, as it prevents the legs to slide backwards from the foot supports, while being tipped.

The leg support is depth adjustable and can be washed in a washing machine.



As an alternative one can choose to get heel stops mounted onto the foot supports.

It prevents the feet to slide backwards of the footplates.

The footplates can still be pushed aside with the heel straps mounted.



HMN since 1934

Technical data

Calf support	
Height on calf support	17 cm
Depth adjustment possible	Yes
Washable at:	80 degrees
Mounted without use of tools	Yes
Made of washable plastic net weave with Velcro tape for adjustment	
Weight	App. 0,3 kg
Article no:	310226

Technical data

Calf support 200 and 300 kg's	
Height on calf support	17 cm
Depth adjustment possible	Yes
Washable at:	80 degrees
Mounted without use of tools	Yes
Made of washable plastic net weave with Velcro tape for adjustment	
Weight	App. 0,3 kg
Article no:	310233

Technical data

Heel straps for foot rests	
Height on heel strap	4,5 cm
Depth adjustment possible	No
Washable at:	80 degrees
Mounted without use of tools	No
Made from washable strap. Nuts and bolts in galvanized steel.	
Weight:	Ca. 0,3 kg
Article no.:	310285



HMN since 1934