



HMN since 1934



PRESSURE-RELIEF SEAT
FOR SHOWER-/COMMODE CHAIRS



DESIGNED TO REDUCE PRESSURE/SHEAR

Pressure-relief seat can treat and prevent pressure ulcers up to grade 3.

HMN's Pressure-relief seat is designed to reduce pressure/shear (displacement) and thereby preventing pressure ulcers. The seat can treat and prevent pressure ulcers up to grade 3.



The anatomical cutout on the seat surface, as well as the flexibility of the material ensures the user stability in the sitting position.

The Pressure-Relief seat only fits on chairs with a horseshoe seat.



The Pressure-Relief seat doesn't fit on M2 300 kg El-Tip, M2 Multi-Tip 250 kg, M2 400 kg, M2 400 kg El-Tip og M2 Mini models.

The seat is placed on the top of the existing seat and can easily be removed.

On the XL models the seat doesn't cover the entire seating surface on the chair.

PRESSURE-RELIEF SEAT

TECHNICAL DATA






If the pillow is cold, it will not have the same softness.




At room temperature, the pillow will have pressure relief properties.

CLEAN WITHOUT ABRASIVES



CHECK THE SEAT FOR DEFECTS AND CRACKS





HMN

since 1934

www.HMN.dk/en

V 1.1 28.03.2025